

—	Event # in original .ht3 file					
—	Event type (0=marker, 1=photon)					
—	Marker or APD input channel #					
—	For marker events:					
	total # of detections since last marker					
	For detection events:					
	Total # of laser pulses emitted thus far					
—	For marker events:					
	sum of time bin values since last marker					
	For detection events:					
	time bin at which detection was received					
1595847	0	01	43689	223482557	(Marker)	
1595848	1	01	153602651	5000	(Detection)	
1595849	1	01	153602665	5067	(Detection)	
1595850	1	01	153602947	5428	(Detection)	
1595852	1	01	153603168	5009	(Detection)	
1595853	1	01	153603230	5064	(Detection)	
⋮						